

# MH-F13 O<sub>2</sub>/DO GAS ANALYZER

DESIGNED FOR PHARMACEUTICAL INDUSTRY



COMPLIES WITH THE REVIEW STANDARDS OF INTERNATIONAL REGULATORY AUTHORITIES

## EQUIPMENT INTRODUCTION

- **Principle:**  
Oxygen molecules quench the fluorescence signal emitted by a specialized sensor spot. The degree of signal attenuation is directly proportional to the oxygen concentration.
- **Methods:**
  1. Invasive: A needle containing the sensor spot pierces the package
  2. Non-Invasive: A fluorescent patch is placed inside a sampling dome, which collects gas from the package

## APPLICATION OBJECTS

- Invasive Needle Method: For vials, infusion bottles, IV Bags, etc
- Non-Invasive Dome Method: For ampoules and other packages unsuitable for puncture

## EQUIPMENT FEATURES

- Suitable for containers with headspace  $\geq 0.1$  mL, including those under negative pressure
- Capable of simultaneous headspace and dissolved oxygen measurement
- Fine-gauge sampling needle allows for easy puncture and is reusable, lowering consumable costs
- Maintenance-free sensor requires no frequent replacement
- Intuitive touchscreen interface
- Unlimited data storage with export/print to PDF or Word
- Software includes database storage, access control, audit trail, and electronic signature functions

## TECHNICAL SPECIFICATIONS

Measurement Range	Headspace O <sub>2</sub> : 0-100% (Full-range sensor)	Dissolved O <sub>2</sub> : 0-44 mg/L (Full-range sensor)
	0-21% (Trace O <sub>2</sub> sensor)	0-9 mg/L (Trace O <sub>2</sub> sensor)
Operating Environment	390 × 350 × 390 mm (LxWxH)	
Accuracy(2-Point Cal.)	Headspace: ±0.2% abs at 20.9%; ±0.02% abs at 1%; ±(2% of reading or 0.01% abs) for trace sensor.	Dissolved: ±0.01 mg/L at 0.44 mg/L; ±0.1 mg/L at 8.8 mg/L
Control	Industrial PC	
Power Supply	100-240 VAC, 50/60 Hz	
Interface	USB	



### Zholion(Shanghai) Technology Co., Ltd.

ADD: 9th Floor, Building11, High-tech Park, Songjiang District, Shanghai  
 WEB: [www.zholionccit.com](http://www.zholionccit.com)  
 EMAIL: [marketing@zholion.com](mailto:marketing@zholion.com)  
 WeChat: 18116456496  
 WhatsApp: +8613331847821